BIOMEK WORKSTATIONS AND ECHO LIQUID HANDLERS IN GENOMIC APPLICATIONS





Echo 525 Genomic applications vary in process complexity and throughput but share the need for high quality and consistent data. Quality is built into the Beckman Coulter Life Sciences instruments to enable your research, whether you need robust automation for a few samples, or need a truly high-throughput solution to increase your capacity and help lower the cost per sample, as well as decreasing the turnaround times.

Biomek Workstations are configurable to handle workflow needs from microliter to milliliter transfers for nucleic acid isolations from primary sample tubes, to preparing sequence-ready libraries. The Echo Acoustic Liquid Handlers (LHs) nanoliter to microliter transfers allow the reduction of reaction costs through reaction volume miniaturization and acoustic non-contact transfers.



1 Sample Prep

Biomek Workstation and Beckman bead-based reagents to automate DNA and RNA extraction.

2 Sample QC

Use either a Biomek Workstation or a combination of Biomek Workstation and Echo LH, depending on your workflow needs.

3 DNA/RNA Normalization

Use either a Biomek Workstation or an Echo LH depending on your workflow needs.

4 Library Preparation

For standard volumes use a Biomek Workstation; or use combination of Echo LH and Biomek Workstation to reduce reaction volumes of precious reagents using the Echo LH, while doing larger volume transfers on the Biomek Workstation.

5 Library Bead Cleanup

Use a Biomek Workstation for 96- and 384-well bead clean ups.

6 Library QC

Use either a Biomek Workstation, or an Echo LH and Biomek Workstation depending on your workflow needs.

Normalization and Pooling

Use a Biomek Workstation for preparing the dilution plates and pooling of the libraries. Alternatively, use an Echo LH for "simultaneous normalization and pooling" for up to 97% faster processing time.

8 Sequencing

Load the pooled libraries onto a flow cell of a sequencer.

Discovery In Motion

Biomek Workstation Solutionsfor Genomic Applications

Biomek 4000 Automated Genomic Workstation



- Extended scalability by 12 deck positions with 4 off-deck mounts
- Flexibility provided by single or 8-channel pipette head options and on the fly automated exchange between tools

Biomek i5 Automated Genomic Workstations



- 25 deck positions to support mid throughput workflows
- Multichannel or Span-8 with gripper options

Biomek i7 Automated Genomic Workstations



- Spacious, open-platform with 45 deck positions for high throughput applications
- Hybrid Workstation with Multichannel pod, Span-8 pod and two independent grippers. Other pipetting configurations are possible as well.

Biomek i-Series Hardware Features



60, 300 and 1200 µL Multichannel heads with 96 and 384 options for flexible transfer options including individual tip(s), column(s), row(s), patterns



Span-8 pod equipped with conductive Liquid Level Sensing (LLS), Allows sample volume tracking through Biomek software



"High-access" rotating gripper optimizes access to high-density decks, allows for direct device integrations.



Optional HEPA filter integration for environmental control. Optional enclosure to provide dust protection.



Integrated devices for process control - e.g. peltiers, shakers and thermocyclers

Successful Library Preparation Evident by QC and Sequencing Data

Obtaining reliable and repeatable results that meet the specifications of reagent vendor standards is crucial to successful automation. Therefore, we partner with customers and leading technology vendors to develop and deliver the best solutions to suit your needs. Our portfolio of Biomek-automated methods are demonstrated to generate quality data using real-world samples.



Biomek i-Series Automated Data Demonstrated Methods

QC Data

Bioanalyzer traces showing consistent sizing between automated and manual libraries. The average size of the manual libraries was 267 bp, while the automated library average size was 266 bp.

E.g. Biomek i7 automated Roche KAPA HyperPrep/HyperPlus method



Sequencing Data

Greater than 97 percent of sequenced reads that mapped to the reference genome.

E.g. Biomek i7 Automated Roche KAPA mRNA HyperPrep Kit



Sequencing Data

Low percentage of ribosomal RNA content in Biomek prepared libraries indicates an efficient capture of mRNA molecules.

E.g. Biomek i7 Automated Roche KAPA mRNA HyperPrep Kit

Complete Your Workflows With Our Genomic Reagents

Biomek automation can be complemented by a growing portfolio of genomic reagent kits from Beckman Coulter Life Sciences, which currently includes nucleic acid extraction and purification solutions for a range of input material. Our high-performance SPRI technology uses paramagnetic beads to selectively immobilize nucleic acids by type and size, and optimized binding conditions enable highly specific separation and cleanup protocols. Portfolio examples include:

Nucleic Acid Purification



AMPure XP

- Removes primers, unincorporated dNTPs, primer dimers, salts and other contaminants
- Eliminates fragments < 50 bp
- Recovers both double- and single-stranded DNA templates
- No PCR degradation after storage at 4° C for 7 days
- Captures DNA > 100 bp

Size Selection



SPRIselect

- Used for fragment size selection for library construction in NGS processes
- Validated for target fragments between 150 bp to 800 bp in length
- Tight manufacturing specifications ensure run to run and lot to lot reproducibility

RNA Isolation



RNAdvance Family

Blood, Cell, Tissue and Viral kits

- Extract RNA from blood (PAXgene tubes), cultured eukaryotic cells, tissue, or saliva and swab transport media
- Produce high-quality RNA compatible with a variety of analysis techniques, such as NGS, microarray, or qRT-PCR

cfDNA Extraction



Apostle MiniMax[™]

High Efficiency cfDNA Isolation Kit

- Isolates cfDNA from 1-5 mL of plasma for liquid biopsy
- Demonstrated compatibility with a variety of collection tubes

Solution Sector Content of Conten

Echo acoustic liquid handling (LH) technology revolutionizes genomic applications by using sound energy to enable highly accurate, fully automated, non-contact dispensing of fluids in nL to µL volumes. By reducing reaction volumes, the Echo LH offers significant cost-savings in reagents, and considerable amount of time-savings for barcoding/indexing and normalization and pooling while preserving data quality.

How Does the Echo LH Move Liquids with Sound?

The heart of the Echo technology is a transducer which receives radio frequency energy and converts it to ultrasounds, which then are transferred to the liquid inside the source well. Using the proprietary dynamic fluid analysis algorithm, the Echo LH determines fluid composition and height, and the power needed to eject a precise volume of fluid into an inverted destination well. This analysis happens in milliseconds, enabling very fast, precise and accurate transfers of nanoliter droplets (less than 8% CV and less than 10% inaccuracy). The desired transfer volumes are achieved by rapidly transferring multiple droplets per second. Transfers can be made from any well in the source to any well in the destination in under a second.





Echo 650 Liquid Handler

Echo 525 Liquid Handler



Acoustic droplet ejection - The transducer is positioned below the source microplate well and emits focused sound energy repeatably to the meniscus of the fluid to be transferred. A stream of 2.5 or 25 nL droplets (model dependent) is reliably ejected into a well of an inverted destination microplate

Low-Cost, Highly Efficient NGS Library Preparation

As the cost of sequencing has gone down, the cost of library preparation is becoming the bottleneck in many high-throughput NGS applications, such as high-throughput amplicon sequencing, or shotgun sequencing in microbial or single cell research. One way to reduce library prep cost is to reduce reaction volumes, therefore there is a need for a liquid handler than can go to sub microliter volumes accurately and precisely. Tip-less, Echo liquid handling systems offer up to 100-fold less reagent costs through miniaturization while ensuring high library quality and repeatability.¹



Integrating Biomek Workstation and Echo Liquid Handler

The use of the Biomek Workstation with an Echo 525 instrument together provide the full workflow automation with the advantages of the acoustic dispensing. The Biomek can directly integrate the Echo LH and control Echo application protocols, allowing the Biomek Workstation and Echo LH to work together to accomplish a complex workflow. This is further enabled by integrating mixing and thermalcycling onto the Biomek to have full workflow capability.

Advantages of an Echo Liquid Handler Integrated to Biomek Workstation

Directly integrate the Echo LH to enable Biomek methods to process plates using Echo Applications version 1.8 protocols.

- Protocols developed using the Echo LH as a stand-alone instrument can be automated using the Biomek Methods
- Run one or more Echo protocols in a Biomek Methods
- Run a full Echo protocol at one time, or run only selected plates within the Echo protocol
- Full use of all features in the Echo application software including identical well content and Echo applications reports
- Visually assign plates to the Echo protocols or use barcoded plates with an Echo picklist transfer file
- Echo LH well volume, dispense and exception data is assigned to plates handled in the Biomek Method for additional calculations and reporting



	BECKMAN			Biomek 15
	 •	0	0	·)
		A 2	10 M	
			1.1	
		•		

	the second s
0	Pun Echo Drotocol
	Run Luio Protocol

elect Echo dev

ECHO v				
C:\Users\Beckman Coulter\Documents\Test2.epr				
Test				
		Plates		
the protocol to run:		Select each plate's location:		
:		Source[1]	P10 ~	
from Source[1] to Destination[1]		Source[2]	P14 ~	
from Source[2] to Destination[1]		Destination[1]	P19 ~	
		Bulk Fills		

Reduced Reaction Volume for Cost and Time Savings



Metagenomic NGS libraries were prepared at 10-times smaller reaction volumes, using the Echo 525 liquid handler. Final libraries produced by the full volume and miniaturized protocols have similar fragment distributions. (A) Bioanalyzer trace of a 5 ng HeLA RNA final library prepared with (A) full volume protocol using manual methods, (B) the miniaturized protocol on Echo 525 LH.²

Mayday MY, Khan LM, Chow ED, Zinter MS, DeRisi JL (2019) Miniaturization and optimization of 384-well compatible RNA sequencing library preparation. PLOS ONE 14(1): e0206194. https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0206194 is licensed under https://creativecommons.org/licenses/by/4.0/ Images used with no changes

Cost savings of reduced reaction volume		FULL VOLUME PROTOCOL	MINIATURIZED PROTOCOL
	Reagents	\$13,427	\$1824
	Other reagents	\$2024	\$282
	Tips	\$1085	\$1020
	Plates/tubes	\$112	\$14
	Source plate	N/A	\$21
	TOTAL	\$16,648	\$3,161

Time savings of reduced reaction volume

5		FULL VOLUME PROTOCOL	MINIATURIZED, AUTOMATED PROTOCOL	
	Calculations/planning	120 minutes	10 minutes	
	Aliquoting RNA	120 minutes	120 minutes	
	Drying RNA	N/A	25 minutes	
	Master mix preparation	120 minutes	10 minutes	
	Dispensing master mixes	215 minutes	45 minutes	
	Incubations	5808 minutes	242 minutes	
	Bead cleanups	2880 minutes	120 minutes	
	Quality check	Using the fragment analyzer and iSeq, QC time is equivalent for 384 samples		
	Pooling	720 minutes	20 minutes	
	TOTAL	9,883 minutes (~166 hours)	592 minutes (~10 hours)	

~80% cost savings and ~90% time savings for reduced reaction volumes (10 times smaller) on the Echo LH compared to the manual and full volume methods.

A General qPCR Workflow



Automated qPCR on Biomek: Flexibility to Automate to Meet Your Throughput and Walkaway Needs

Automation on Biomek Workstations replaces repetitive pipetting steps in PCR/qPCR setup. Biomek allows the flexibility to work with multiple plate formats including 96 and 384, with any combination of reaction components. On-deck temperature control by Peltier maintains sensitive reagents at the appropriate temperature, minimizing the possibility of assay failure. Complete workflow automation is achieved through integration of automated thermocyclers, reducing hands-on time and associated human errors.

Biomek i7 Automated GAPDH qPCR Assay

The calculated CT values are in excellent agreement between manual reaction preparation (21.60 \pm 0.06) and Biomek automated reaction prep (21.68 \pm 0.08) for GAPDH qPCR assay. As expected, a concentration dependent shift in the fluorescence trace was observed. The highest amount of RNA input (10 ng) produced a CT value of 18.2 \pm 0.1, while the lowest concentration tested (3 pg) had a CT value of 30.33 \pm 0.09. When the calculated CT value was plotted as a function of RNA input an excellent linearity (R2 = 0.9986) was observed.



Integrating the Echo LH into an Access System for Robotic Plate Handling



An example Access Workstation with third-party devices for peeling, centrifugation, bulkfilling, sealing and gPCR cycling, to automate reaction setup for qPCR screening around an Echo Liquid Handler.

The benefits of using the Echo Liquid Handler can be a solution for a number of assays. The Access Systems are designed to be modular and configurable workflow automation. Using robotic plate handling and third party devices, the systems optimize the workflow scheduling around the Echo protocols. Integrate two Echo LHs to allow simultaneous use of different instruments to assemble various assays. The Access Workstation is a benchtop robotic platform for most any application. And the Access Single Robot System and Dual Robot System (Access SRS and DRS) are designed with docking carts where many third party devices are integrated vertically in a small footprint for greater capability, with optional HEPA air handling or low humidity HVAC attachment for the appropriate environment in the enclosure for your assays.

RT-qPCR on Echo 525 LH: 0.5 μ L RT-qPCR, Similar Performance to 12 μ L

Effect of reference compound A on hepatitis C virus (HCV) replicon RNA synthesis in human hepatoma (Huh7) cells was studied using RTqPCR (TaqMan) analysis of the HCV replicon RNA level (C_{τ}). Huh7 cells were treated with increasing concentrations of reference compound A for 72 h. Increases in C_{T} values indicate decreasing replicon RNA levels; each C_{τ} reflects a two-fold change in RNA level from baseline. Same performance was observed between reduced volume reactions, as compared to standard reaction volumes, set up using traditional methods.

MM: master mix with primers and probe CL: cells with lysis buffer



35

Agrawal S, Cifelli S, Johnstone R, Pechter D, Barbey DA, Lin K, Allison T, Agrawal S, Rivera-Gines A, Milligan JA, Schneeweis J, Houle K, Struck AJ, Visconti R, Sills M, Wildey MJ, Utilizing Low-Volume Aqueous Acoustic Transfer with the Echo 525 to Enable Miniaturization of qRT-PCR (Vol. 21(1)) pp 57-63, copyright © 2015 by Society for Laboratory Automation and Screening, Reprinted by Permission of SAGE Publications

Integrating the Biomek Workstation into a Complete, Multi-Device Robotic System

Genomic applications using the Biomek Workstation can include robotic plate handling to obtain full walkaway solutions, for sample to results processing and minimizing errors that can come from manual processing. The engineers in the Access Solution Group at Beckman Coulter Life Sciences can help create a solution to meet these needs, with the ability to integrate more than 300 third-party devices from over 60 manufacturers in different laboratory settings. Our SAMI EX scheduling software works to plan workflows and batching to maintain incubation times for reproducibility of results.



An example laboratory system with Biomek i7 Workstation and integrated devices to fully automate nucleic acid isolation from sample tubes, amplification reaction setup and gPCR/RTgPCR automation, in one system. Users load the sample tubes and needed consumables, and the system processes the samples to generate results.

Consumables Biomek and Echo Qualified for Genomic Applications





Biomek Pipetting Tips

- **RNase- and DNase-free*:** Ensures high-quality nucleic acid purification and reliable testing results.
- DNA-free (human & mouse) and PCR inhibitor-free: Confirms the absence of of any contaminants (microbial or mammalian) that contribute to erroneous results and interfere with PCR.
- **Pyrogen-free:** Addresses pyrogen-sensitive applications, such as cell transformations and preparative protocols that use extracted biological material for in vivo experimentation.
- **Trace metal-free:** Reduces chemical interference due to chelation and denaturation.

* Free means: Bio-certification "free of" claims are defined as the lower limit of detection based on the sensitivity of the test method or instrumentation used.

Echo Qualified Source Microplates ensure precise, accurate liquid transfers, delivering the best possible results in miniaturized assays. Echo qualified source microplates are deionized and are packaged in anti-static bags to ensure precise drop placement. Microplates are available with low dead volume, custom barcode, DNAse and RNAse free or sterile options.





The MicroClime Environmental Lid minimizes edge effects and preserves the concentrations of solutions in microplate wells. The novel fluid-absorbing matrix inside the lid can be loaded with water or a compatible solvent to create a vapor barrier, protecting the sample from evaporation and greatly reducing edge effects, and therefore increasing the microplate efficiency by 37%.



The Brooks Life Sciences FluidX AcoustiX Sample Tubes transform workflows by enabling acoustic dispensing directly from tubes, on Echo 650 Series LHs. The tubes can be accessed individually – ideal for applications that require subsets of large reagent or sample libraries to be accessed frequently.

Beckman Coulter Life Sciences Service and Support



Maximize Your Instrument Performance

Timely service and preventive maintenance are essential for optimal instrument performance and data quality. Beckman Coulter Life Sciences offers a wide range of service contract options for the Biomek Workstations, Echo Liquid Handlers and Access Laboratory Workstations and Systems to fit every lab's needs and budget.

Our application scientists have extensive experience with genomic workflows and provide support and training on all our product offerings to help you maximize the use of your systems.

beckman.com

REFERENCES

- Shapland EB, Holmes V, Reeves CD, et al. Low-Cost, High-throughput Sequencing of DNA Assemblies Using a Highly Multiplexed Nextera Process. Shapland et al., ACS Synth. Biol., 2015. doi.org/10.1021/sb500362n
- Mayday MY, Khan LM, Chow ED, Zinter MS, DeRisi JL (2019) Miniaturization and optimization of 384-well compatible RNA sequencing library preparation. PLoS ONE 14(1): e0206194. https://doi.org/10.1371/journal.pone.0206194

The method/application illustrated is for demonstration only, and is not validated by Beckman Coulter. Beckman Coulter makes no warranties of any kind whatsoever express or implied, with respect to this protocol, including but not limited to warranties of fitness for a particular purpose or merchantability or that the protocol is non-infringing. All warranties are expressly disclaimed. Your use of the method is solely at your own risk, without recourse to Beckman Coulter. Not intended or validated for use in the diagnosis of disease or other conditions.

Products identified are not for use in diagnostics procedures.

© 2021 Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries. ECHO is a trademark or registered trademark of Labcyte Inc. in the United States and other countries.



All other product names and brands are properties of their respective owners.

For Beckman Coulter's worldwide office locations and phone numbers, please visit Contact Us at beckman.com 21.08.2274.AUTO